

**Sterile Saline Eye Wash Solution**

**Sterile Eyewash**

**1 Identification**

**1.1 Product Code Description**

2142	20ml Blue Dot Sterile Emergency Eye Wash Pods (Each)
2144	250ml Blue Dot Emergency Eye Wash Solution (Each)
2143	500ml Blue Dot Eye Wash Solution in round bottle (each)
2403008	500ml Astroplast Sterile Eye Wash Solution (Each)
1047009	500ml Astroplast Sterile Eye Wash Solution (case 3 x Box 3) (11663WC)
1047091	500ml Astroplast Sterile Eye Wash Solution (case 10) (11656WC)
2401020	500ml Astroplast Sterile Eye Wash Solution (Box 3)
2404091	500ml Astroplast Sterile Eye Wash Solution (Case 10) Innovative Safety
2404085	500ml Astroplast Sterile Eye Wash Solution (Case 3 x Box 2) (10848) (347490)
2404034A	500ml Astroplast Sterile Eye Wash Solution VOSS (Case 10) Export
2401023	500ml Astroplast Twist & Open Large Eye Wash Station Complete (Each)
2403011	500ml Astroplast Twist & Open Eye Wash Station Complete (Each)
2404034	500ml Twist & Open Eye Wash Station Complete (Each)
2404086	500ml Astroplast Twist & Open Eye Wash Solution Complete (Box 3) UK
2404087	500ml Astroplast Twist & Open Sterile Eye Wash Solution (Box 3) EUR
2404088	500ml Astroplast Twist & Open Sterile Eye Wash Solution with Euro Cap (Case 10)
2404094	500ml Astroplast Twist & Open Sterile Eye Wash Solution (Case 3 x Box 3) (6608875) Lyreco
2404093	500ml Astroplast Twist & Open Sterile Eye Wash Solution (Case 4 x Box 2) (13346) (543408)
2405105	500ml Astroplast Twist & Open Sterile Solution New Round PK3 +2 x Free Skate (Each)
2143SP	500ml SP Services Eye Wash Solution in Round Bottle (Each)
2404089	Astroplast Innovative Safety 500ml Eye Wash (English & French)

## Sterile Saline Eye Wash Solution

### 1.2 Identified use of the substance

- Eye wash is recommended for rinsing foreign bodies (such as dust, dirt, metal and wood splinters) and (chemicals such as solvents and oil) from the eyes.
- Position container above the eye and gently squeeze so that the fluid flows freely across the eye towards the outer corner.

### 2 Classification

- Class 1 Sterile, according to Rule 5 of Annex IX, MDD93/42/EEC
- Not a hazardous substance or mixture according to regulation (EC)No. 1272/2008

### 3 Sterilisation

- Filtration

### 4 Information on ingredients

Component	%	CAS #	EC #	REACH Reg. No.
Sodium Chloride	0.9	7647-14-5	231-598-3	N/A

### 5 First Aid Measures

N/A

#### 5.1 Indication of any immediate Medical attention and Special treatment

N/A

### 6 Extinguishing Media

Only the packaging may cause a fire, in the event of a large fire breathing apparatus to be worn.

- Suitable extinguishing media: Foam. Dry powder. Carbon dioxide. Water spray. Sand.
- Unsuitable extinguishing media: Do not use a heavy water stream.

### 7 Accidental Release

N/A

### 8 Handling and Storage

- Keep out of children's reach
- To be stored between 5°C- 35°C

### 9 Appearance

- Sterile Saline irrigation is a clear colourless solution

**Sterile Saline Eye Wash Solution**

## 10 Stability

- Product is stable

## 11 Disposal

- Saline solution can be poured into any drain

## 12 Transportation

- The product is not Hazardous for transport purposes